



## Membrane pump M401

P/N YYM401

is a pump for gasses pumping to do underpressure and overpressure for lab use.


The parts of pump where are in contact with pumped gas are manufacture from chemical resist materials as stainless steel, PTFE and hipalon so that it can be used for aggressive gasses pumping.

### Sample of use:

- by mobile phase filtration using SolVac Filter Holder
- by Solid Phase Extraction

### Specification

Pumped volume at zero pressure	6 NI /min ± 10%
Pumped volume at pressure of 50kPa	3,5 NI /min ± 10%
Pumped volume at pressure of 100kPa	1,8 NI /min ± 10%
Nominal voltage	230 V ± 10% , 50 Hz
Nominal input	65 VA
Dimensions (W x H x D)	98 x 140 x 145 mm (3,86" x 5,51" x 5,70)
Weight	3,1 kg



Informative power curve of pump

☺ underpressure

overpressure ☹